



Ankle-Arm Blood Pressure

Id#: IDNO

Acrostic: \_\_\_\_\_

Date:   /   /      
 Month Day Year

**Pulse Obliteration Pressure:**

*Inflate the cuff slowly until the pulse is no longer audible. Inflate to 30 mm Hg above the level at which pulse sound disappeared.*

**1 Systolic Readings:**

Record in this order

	Systolic
Right brachial	<input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <b>RBRACH1</b>
Right dorsalis pedis	<input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <b>RDPEDIS1</b>
Right posterior tibial	<input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <b>RPTIB1</b>
	Systolic
Left dorsalis pedis	<input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <b>LDPEDIS1</b>
Left posterior tibial	<input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <b>LPTIB1</b>
Left brachial	<input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <b>LBRACH1</b>

- 2 Procedure was:** **AABPCMP1**
- Completed successfully → Skip to Technician ID
  - Not completed → Continue to #3

**3 Reason procedure was not completed:**

<b>RIGHT SIDE:</b>	<b>YES</b>	<b>NO</b>
Unable to occlude	<input type="radio"/>	<input type="radio"/> <b>RLUOCCL1</b>
Ulceration	<input type="radio"/>	<input type="radio"/>
Amputation	<input type="radio"/>	<input type="radio"/>
Unable to locate DP	<input type="radio"/>	<input type="radio"/> <b>RULOCDP1</b>
Unable to locate PT	<input type="radio"/>	<input type="radio"/> <b>RULOCPT1</b>

Other:

<b>LEFTSIDE:</b>	<b>YES</b>	<b>NO</b>
Unable to occlude	<input type="radio"/>	<input type="radio"/> <b>LLUOCCL1</b>
Ulceration	<input type="radio"/>	<input type="radio"/>
Amputation	<input type="radio"/>	<input type="radio"/>
Unable to locate DP	<input type="radio"/>	<input type="radio"/> <b>LULOCDP1</b>
Unable to locate PT	<input type="radio"/>	<input type="radio"/> <b>LULOCPT1</b>

Other:

Technician ID#:

Reviewer ID#:

Data Entry ID #: